

**Ptma Rabbit anti-Human Polyclonal (aa8-39) (PE) Antibody - LS-C267527 - LSBio**

<b>CatalogID:</b>	LS-C267527
<b>Target:</b>	prothymosin, alpha (Ptma)
<b>Synonyms:</b>	Ptma Antibody, Gene sequence 28 Antibody, Prothymosin, alpha Antibody, TMSA Antibody, Prothymosin alpha Antibody, Prothymosin alpha protein Antibody, PTMA Antibody
<b>Host</b>	Ptma antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	Ptma antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	Ptma antibody was raised against kLH-conjugated synthetic peptide (aa8-39) from the N-terminal region of human PTMA.
<b>Specificity:</b>	Human PTMA
<b>Epitope:</b>	aa8-39
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences