

PTPRCAP / CD45 AP Rabbit anti-Human Polyclonal (aa89-118) (AP) Antibody - LS-C267570 - LSBio

CatalogID:	LS-C267570
Target:	protein tyrosine phosphatase, receptor type, C-associated protein (PTPRCAP)
Synonyms:	PTPRCAP Antibody, CD45-associated protein Antibody, LPAP Antibody, CD45 associated protein Antibody, CD45-AP Antibody, PTPRC-associated protein Antibody
Host	PTPRCAP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	PTPRCAP / CD45 AP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PTPRCAP / CD45 AP antibody was raised against kLH-conjugated synthetic peptide (aa89-118) from internal region of human PTPRCAP.
Specificity:	Human PTPRCAP / CD45 AP.
Epitope:	aa89-118
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences