

**C6orf108 Rabbit anti-Human Polyclonal (aa27-56) (FITC) Antibody - LS-C267999 - LSBio**

<b>CatalogID:</b>	LS-C267999
<b>Target:</b>	2'-deoxynucleoside 5'-phosphate N-hydrolase 1 (DNPH1)
<b>Synonyms:</b>	DNPH1 Antibody, C6orf108 Antibody, DJ330M21.3 Antibody, RCL Antibody, RP3-330M21.3 Antibody, C-Myc-responsive protein Rcl Antibody, Putative c-Myc-responsive Antibody
<b>Host</b>	DNPH1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	C6orf108 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	C6orf108 antibody was raised against kLH-conjugated synthetic peptide (aa27-56) from the N-terminal region of human RCL.
<b>Specificity:</b>	Human C6orf108.
<b>Epitope:</b>	aa27-56
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences