

C6orf108 Rabbit anti-Human Polyclonal (aa27-56) (Azide-free) (HRP) Antibody - LS-C268000 - LSBio

CatalogID:	LS-C268000
Target:	2'-deoxynucleoside 5'-phosphate N-hydrolase 1 (DNPH1)
Synonyms:	DNPH1 Antibody, C6orf108 Antibody, DJ330M21.3 Antibody, RCL Antibody, RP3-330M21.3 Antibody, C-Myc-responsive protein Rcl Antibody, Putative c-Myc-responsive Antibody
Host	DNPH1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	C6orf108 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	C6orf108 antibody was raised against kLH-conjugated synthetic peptide (aa27-56) from the N-terminal region of human RCL.
Specificity:	Human C6orf108.
Epitope:	aa27-56
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, azide-free.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences