

RELA / NFkB p65 Rabbit anti-Human Polyclonal (pSer536) (FITC) Antibody - LS-C268076 - LSBio	
CatalogID:	LS-C268076
Target:	v-rel avian reticuloendotheliosis viral oncogene homolog A (RELA)
Synonyms:	RELA Antibody, NF-kappa-B p65delta3 Antibody, NFkB3 Antibody, p65 Antibody, Transcription factor p65 Antibody
Host	RELA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	RELA / NFkB p65 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	RELA / NFkB p65 antibody was raised against KLH-conjugated synthetic phosphopeptide around Ser536 of human RELA.
Specificity:	Human RELA / NFkB p65
Epitope:	pSer536
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	