

RHBG Rabbit anti-Human Polyclonal (aa421-450) (PE) Antibody - LS-C268216 - LSBio

CatalogID:	LS-C268216
Target:	Rh family, B glycoprotein (gene/pseudogene) (RHBG)
Synonyms:	RHBG Antibody, Ammonium transporter Rh type B Antibody, Rh family type B glycoprotein Antibody, SLC42A2 Antibody, Rh type B glycoprotein Antibody
Host	RHBG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RHBG antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RHBG antibody was raised against kLH-conjugated synthetic peptide (aa421-450) from the C-terminus of human RHBG.
Specificity:	Human RHBG.
Epitope:	aa421-450
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences