

RIPK3 / RIP3 Rabbit anti-Human Polyclonal (aa479-509) (FITC) Antibody - LS-C268284 - LSBio	
CatalogID:	LS-C268284
Target:	receptor-interacting serine-threonine kinase 3 (RIPK3)
Synonyms:	RIPK3 Antibody, Receptor interacting protein 3 Antibody, RIP-like protein kinase 3 Antibody, RIP3 Antibody, RIP-3 Antibody, Receptor-interacting protein 3 Antibody
Family / Subfamily:	Protein Kinase / RIP
Host	RIPK3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	RIPK3 / RIP3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RIPK3 / RIP3 antibody was raised against kLH-conjugated synthetic peptide (aa479-509) from the C-terminus of human RIPK3.
Specificity:	Human RIPK3 / RIP3
Epitope:	aa479-509
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	