

RNF185 Rabbit anti-Human Polyclonal (aa91-121) (FITC) Antibody - LS-C268330 - LSBio

CatalogID:	LS-C268330
Target:	ring finger protein 185 (RNF185)
Synonyms:	RNF185 Antibody, BSK65-MONO1 Antibody, BSK65-PANC1 Antibody, BSK65-PANC2 Antibody, BSK65-TEST1 Antibody, BSK65-TEST2 Antibody, BSK65-TEST3 Antibody, Ring finger protein 185 Antibody, BSK65-MONO2 Antibody
Host	RNF185 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	RNF185 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RNF185 antibody was raised against kLH-conjugated synthetic peptide (aa91-121) from internal region of human RN185.
Specificity:	Human RNF185
Epitope:	aa91-121
Reactivity:	Human, Mouse, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences