

RNASEK Rabbit anti-Human Polyclonal (aa34-64) (PE) Antibody - LS-C268349 - LSBio

CatalogID:	LS-C268349
Target:	ribonuclease, RNase K (RNASEK)
Synonyms:	RNASEK Antibody, Ribonuclease kappa Antibody, RNase K Antibody, Ribonuclease, RNase K Antibody, RNase kappa Antibody
Host	RNASEK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RNASEK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RNASEK antibody was raised against kLH-conjugated synthetic peptide (aa34-64) from internal region of human RNASEK.
Specificity:	Human RNASEK
Epitope:	aa34-64
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences