

RYP Rabbit anti-Human Polyclonal (aa168-198) (PE) Antibody - LS-C268690 - LSBio

CatalogID:	LS-C268690
Target:	receptor-like tyrosine kinase (RYK)
Synonyms:	RYK Antibody, D3S3195 Antibody, Hydroxyaryl-protein kinase Antibody, JTK5 Antibody, JTK5A protein tyrosine kinase Antibody, RYK1 Antibody, Vik Antibody, Tyrosine-protein kinase RYK Antibody, Receptor-like tyrosine kinase Antibody, JTK5A Antibody, Tyrosine kinase ryk Antibody
Family / Subfamily:	Protein Kinase / RYK
Host	RYK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RYK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RYK antibody was raised against kLH-conjugated synthetic peptide (aa168-198) within aa 150~200 of human RYK.
Specificity:	Human RYK.
Epitope:	aa168-198
Reactivity:	Human, Mouse, Rat
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences