

**RYK Rabbit anti-Human Polyclonal (aa569-599) (PE) Antibody - LS-C268695 - LSBio**

|                             |   |
|-----------------------------|---|
| <b>CatalogID:</b>           | LS-C268695  |
| <b>Target:</b>              | receptor-like tyrosine kinase (RYK)   |
| <b>Synonyms:</b>            | RYK Antibody, D3S3195 Antibody, Hydroxyaryl-protein kinase Antibody, JTK5 Antibody, JTK5A protein tyrosine kinase Antibody, RYK1 Antibody, Vik Antibody, Tyrosine-protein kinase RYK Antibody, Receptor-like tyrosine kinase Antibody, JTK5A Antibody, Tyrosine kinase ryk Antibody |
| <b>Family / Subfamily:</b>  | Protein Kinase / RYK  |
| <b>Host</b>                 | RYK antibody was produced in Rabbit   |
| <b>Clonality:</b>           | Polyclonal  |
| <b>Isotype:</b>             | IgG   |
| <b>Conjugations:</b>        | Phycoerythrin (PE)  |
| <b>Immunogen Species:</b>   | RYK antibody was raised against Human   |
| <b>Antigen Type:</b>        | Synthetic peptide   |
| <b>Immunogen:</b>           | RYK antibody was raised against kLH-conjugated synthetic peptide (aa569-599) from the C-terminus of human RYK.  |
| <b>Specificity:</b>         | Human RYK.  |
| <b>Epitope:</b>             | aa569-599   |
| <b>Reactivity:</b>          | Human   |
| <b>Purification:</b>        | Affinity purified   |
| <b>Presentation:</b>        | PBS, pH 7.2, 0.09% sodium azide.  |
| <b>Recommended Storage:</b> | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.   |
| <b>Usage Summary:</b>       | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.  |
| <b>Uses:</b>                | IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)  |
| <b>Size:</b>                | 200 µl  |
| <b>Requested From:</b>      | Japan   |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences