

SHBG Rabbit anti-Human Polyclonal (aa68-96) (FITC) Antibody - LS-C269326 - LSBio

CatalogID:	LS-C269326
Target:	sex hormone-binding globulin (SHBG)
Synonyms:	SHBG Antibody, ABP Antibody, SBP Antibody, Sex hormone-binding globulin Antibody, Sex steroid-binding protein Antibody, TEBG Antibody
Host	SHBG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SHBG antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SHBG antibody was raised against KLH-conjugated synthetic peptide (aa68-96) from the N-terminal region of human SHBG.
Specificity:	Human SHBG.
Epitope:	aa68-96
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences