

SOX2 Rabbit anti-Human Polyclonal (pSer246) (AP) Antibody - LS-C270009 - LSBio

| | |
|-----------------------------|--|
| CatalogID: | LS-C270009 |
| Target: | SRY (sex determining region Y)-box 2 (SOX2) |
| Synonyms: | SOX2 Antibody, MCOPS3 Antibody, Transcription factor SOX2 Antibody, ANOP3 Antibody, Transcription factor SOX-2 Antibody, SRY-related HMG-box gene 2 Antibody |
| Family / Subfamily: | Transcription factor |
| Host | SOX2 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Alkaline Phosphatase (AP) |
| Immunogen Species: | SOX2 antibody was raised against Human |
| Antigen Type: | Synthetic phospho-peptide |
| Immunogen: | SOX2 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser246 of human Sox2. |
| Specificity: | Human SOX2. |
| Epitope: | pSer246 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | ICC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences