

CD43 Rabbit anti-Human Polyclonal (aa31-60) (FITC) Antibody - LS-C270237 - LSBio

CatalogID:	LS-C270237
Target:	sialophorin (SPN)
Synonyms:	SPN Antibody, GPL115 Antibody, Leukosialin Antibody, Leukocyte sialoglycoprotein Antibody, Galactoglycoprotein Antibody, GALGP Antibody, Sialophorin Antibody, LSN Antibody, CD43 Antibody, CD43 antigen Antibody
Host	SPN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD43 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD43 antibody was raised against kLH-conjugated synthetic peptide (aa31-60) from the N-terminal region of human SPN/CD43.
Specificity:	Human CD43.
Epitope:	aa31-60
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences