

**S. aureus Delta Hemolysin Rabbit anti-Staphylococcus aureus Polyclonal (FITC) Antibody - LS-C270487 - LSBio**

<b>CatalogID:</b>	LS-C270487
<b>Target:</b>	S. aureus Delta Hemolysin
<b>Host</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Staphylococcus aureus
<b>Antigen Type:</b>	Hapten
<b>Immunogen:</b>	Delta hemolysin protein of Staphylococcus aureus, conjugated to KLH.
<b>Specificity:</b>	Staphylococcus aureus S. aureus Delta Hemolysin.
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences