

TBCA Rabbit anti-Mouse Polyclonal (aa37-66) (PE) Antibody - LS-C271103 - LSBio

CatalogID:	LS-C271103
Target:	tubulin folding cofactor A (TBCA)
Synonyms:	TBCA Antibody, Chaperonin cofactor A Antibody, CFA Antibody, Co-chaperonin 'cofactor a' Antibody, Tubulin-specific chaperone A Antibody, Tubulin-folding cofactor A Antibody, TCP1-chaperonin cofactor A Antibody, Tubulin folding cofactor A Antibody
Host	TBCA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TBCA antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	TBCA antibody was raised against kLH-conjugated synthetic peptide (aa37-66) from internal region of human TBCA.
Specificity:	Mouse TBCA.
Epitope:	aa37-66
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences