

TEC Rabbit anti-Human Polyclonal (aa25-55) (FITC) Antibody - LS-C271254 - LSBio

CatalogID:	LS-C271254
Target:	tec protein tyrosine kinase
Synonyms:	TEC Antibody, Tyrosine-protein kinase Tec Antibody, PSCTK4 Antibody, Tec family Antibody, Tec protein tyrosine kinase Antibody
Family / Subfamily:	Protein Kinase / BTK/TEC
Host	TEC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TEC antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TEC antibody was raised against kLH-conjugated synthetic peptide (aa25-55) from the N-terminal region of human TEC.
Specificity:	Human TEC
Epitope:	aa25-55
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences