

CD30 Rabbit anti-Human Polyclonal (aa470-500) (FITC) Antibody - LS-C271883 - LSBio

CatalogID:	LS-C271883
Target:	tumor necrosis factor receptor superfamily, member 8 (TNFRSF8)
Synonyms:	TNFRSF8 Antibody, CD30L receptor Antibody, CD30 Antibody, Cytokine receptor CD30 Antibody, D1S166E Antibody, Ki-1 Antibody, Ki-1 antigen Antibody, CD30 antigen Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	TNFRSF8 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD30 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD30 antibody was raised against kLH-conjugated synthetic peptide (aa470-500) from the Y479 region of human TNFRSF8.
Specificity:	Human CD30.
Epitope:	aa470-500
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences