

Hamartin / TSC1 Rabbit anti-Rat Polyclonal (pSer1141) (FITC) Antibody - LS-C272333 - LSBio

CatalogID:	LS-C272333
Target:	tuberous sclerosis 1 (TSC1)
Synonyms:	TSC1 Antibody, LAM Antibody, Tuberous sclerosis 1 protein Antibody, Tuberous sclerosis 1 Antibody, Hamartin Antibody, KIAA0243 Antibody, TSC Antibody, Tumor suppressor Antibody
Host	TSC1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Hamartin / TSC1 antibody was raised against Rat
Antigen Type:	Synthetic phospho-peptide
Immunogen:	Hamartin / TSC1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser1141 of rat TSC1.
Specificity:	Rat Hamartin / TSC1
Epitope:	pSer1141
Reactivity:	Rat, Rabbit
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences