

| Hamartin / TSC1 Rabbit anti-Human Polyclonal (pSer682) (PE) Antibody - LS-C272401 - LSBio | |
|---|--|
| CatalogID: | LS-C272401 |
| Target: | tuberous sclerosis 1 (TSC1) |
| Synonyms: | TSC1 Antibody, LAM Antibody, Tuberous sclerosis 1 protein Antibody, Tuberous sclerosis 1 Antibody, Hamartin Antibody, KIAA0243 Antibody, TSC Antibody, Tumor suppressor Antibody |
| Host | TSC1 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | Hamartin / TSC1 antibody was raised against Human |
| Antigen Type: | Synthetic phospho-peptide |
| Immunogen: | Hamartin / TSC1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser682 of human TSC1. |
| Specificity: | Human Hamartin / TSC1 |
| Epitope: | pSer682 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | ELISA, Dot Blot (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/30/2014 | |
| © 2014 LifeSpan BioSciences | |