

Hamartin / TSC1 Rabbit anti-Mouse Polyclonal (pThr1043) (FITC) Antibody - LS-C272405 - LSBio	
CatalogID:	LS-C272405
Target:	tuberous sclerosis 1 (TSC1)
Synonyms:	TSC1 Antibody, LAM Antibody, Tuberous sclerosis 1 protein Antibody, Tuberous sclerosis 1 Antibody, Hamartin Antibody, KIAA0243 Antibody, TSC Antibody, Tumor suppressor Antibody
Host	TSC1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Hamartin / TSC1 antibody was raised against Mouse
Antigen Type:	Synthetic phospho-peptide
Immunogen:	Hamartin / TSC1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Thr1043 of mouse TSC1.
Specificity:	Mouse Hamartin / TSC1
Epitope:	pThr1043
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	