

**Hamartin / TSC1 Rabbit anti-Mouse Polyclonal (pTyr304) (PE) Antibody - LS-C272461 - LSBio**

<b>CatalogID:</b>	LS-C272461
<b>Target:</b>	tuberous sclerosis 1 (TSC1)
<b>Synonyms:</b>	TSC1 Antibody, LAM Antibody, Tuberous sclerosis 1 protein Antibody, Tuberous sclerosis 1 Antibody, Hamartin Antibody, KIAA0243 Antibody, TSC Antibody, Tumor suppressor Antibody
<b>Host</b>	TSC1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	Hamartin / TSC1 antibody was raised against Mouse
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	Hamartin / TSC1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr304 of mouse TSC1.
<b>Specificity:</b>	Mouse Hamartin / TSC1
<b>Epitope:</b>	pTyr304
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences