

Tuberin / TSC2 Rabbit anti-Mouse Polyclonal (pSer1365) (FITC) Antibody - LS-C272524 - LSBio	
CatalogID:	LS-C272524
Target:	tuberous sclerosis 2 (TSC2)
Synonyms:	TSC2 Antibody, LAM Antibody, Tuberous sclerosis 2 protein Antibody, TSC4 Antibody, Tuberin Antibody, Tuberous sclerosis 2 Antibody
Host	TSC2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Tuberin / TSC2 antibody was raised against Mouse
Antigen Type:	Synthetic phospho-peptide
Immunogen:	Tuberin / TSC2 antibody was raised against KLH-conjugated synthetic phosphopeptide around Ser1365 of mouse TSC2.
Specificity:	Mouse Tuberin / TSC2
Epitope:	pSer1365
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	