

Tuberin / TSC2 Rabbit anti-Rat Polyclonal (pThr1373) (FITC) Antibody - LS-C272556 - LSBio	
<b>CatalogID:</b>	LS-C272556
<b>Target:</b>	tuberous sclerosis 2 (TSC2)
<b>Synonyms:</b>	TSC2 Antibody, LAM Antibody, Tuberous sclerosis 2 protein Antibody, TSC4 Antibody, Tuberin Antibody, Tuberous sclerosis 2 Antibody
<b>Host</b>	TSC2 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Tuberin / TSC2 antibody was raised against Rat
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	Tuberin / TSC2 antibody was raised against KLH-conjugated synthetic phosphopeptide around Thr1373 of rat TSC2.
<b>Specificity:</b>	Rat Tuberin / TSC2
<b>Epitope:</b>	pThr1373
<b>Reactivity:</b>	Rat, Rabbit
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/30/2014 © 2014 LifeSpan BioSciences	