

UBA1 / UBE1 Rabbit anti-Human Polyclonal (aa1-30) (FITC) Antibody - LS-C272828 - LSBio

CatalogID:	LS-C272828
Target:	ubiquitin-like modifier activating enzyme 1 (UBA1)
Synonyms:	UBA1 Antibody, A1S9 Antibody, A1ST Antibody, A1S9T Antibody, AMCX1 Antibody, GXP1 Antibody, Protein A1S9 Antibody, SMAX2 Antibody, Ubiquitin-activating enzyme E1 Antibody, UBE1X Antibody, POC20 Antibody, UBE1 Antibody, UBA1A Antibody
Host	UBA1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	UBA1 / UBE1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	UBA1 / UBE1 antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminal region of human UBA1.
Specificity:	Human UBA1 / UBE1
Epitope:	aa1-30
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 8/28/2014

© 2014 LifeSpan BioSciences