

UPRT Rabbit anti-Mouse Polyclonal (aa203-233) (FITC) Antibody - LS-C273085 - LSBio

CatalogID:	LS-C273085
Target:	uracil phosphoribosyltransferase (FUR1) homolog (<i>S. cerevisiae</i>) (UPRT)
Synonyms:	UPRT Antibody, FUR1 Antibody, UMP pyrophosphorylase Antibody, UPP Antibody, RP11-311P8.3 Antibody, UPRTase Antibody
Host	UPRT antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	UPRT antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	UPRT antibody was raised against kLH-conjugated synthetic peptide (aa203-233) from the C-terminus of human UPRT.
Specificity:	Mouse UPRT.
Epitope:	aa203-233
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences