

VEGFB Rabbit anti-Human Polyclonal (aa117-147) (FITC) Antibody - LS-C273246 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C273246 |
| Target: | vascular endothelial growth factor B (VEGFB) |
| Synonyms: | VEGFB Antibody, VRF Antibody, VEGF-B Antibody, VEGF-related factor Antibody, VEGFL Antibody |
| Family / Subfamily: | PDGF / not assigned-PDGF |
| Host | VEGFB antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | VEGFB antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | VEGFB antibody was raised against kLH-conjugated synthetic peptide (aa117-147) from internal region of human VEGF2. |
| Specificity: | Human VEGFB. |
| Epitope: | aa117-147 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences