

VEGFB Rabbit anti-Human Polyclonal (aa117-147) (PE) Antibody - LS-C273248 - LSBio

CatalogID:	LS-C273248
Target:	vascular endothelial growth factor B (VEGFB)
Synonyms:	VEGFB Antibody, VRF Antibody, VEGF-B Antibody, VEGF-related factor Antibody, VEGFL Antibody
Family / Subfamily:	PDGF / not assigned-PDGF
Host	VEGFB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	VEGFB antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	VEGFB antibody was raised against kLH-conjugated synthetic peptide (aa117-147) from internal region of human VEGF2.
Specificity:	Human VEGFB.
Epitope:	aa117-147
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences