

XKRX Rabbit anti-Human Polyclonal (aa316-345) (FITC) Antibody - LS-C273714 - LSBio

CatalogID:	LS-C273714
Target:	XK, Kell blood group complex subunit-related, X-linked (XKRX)
Synonyms:	XKRX Antibody, Membrane protein XPLAC Antibody, XK-related protein 2 Antibody, XRG2 Antibody, XKR2 Antibody, XPLAC Antibody
Host	XKRX antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	XKRX antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	XKRX antibody was raised against kLH-conjugated synthetic peptide (aa316-345) from the C-terminus of human XKRX.
Specificity:	Human XKRX
Epitope:	aa316-345
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences