

YAP / YAP1 Rabbit anti-Human Polyclonal (Ser127) (PE) Antibody - LS-C273773 - LSBio

CatalogID:	LS-C273773
Target:	Yes-associated protein 1 (YAP1)
Synonyms:	YAP1 Antibody, 65 kDa Yes-associated protein Antibody, YAP Antibody, Yes-associated protein beta Antibody, Yes-associated protein delta Antibody, Yorkie homolog Antibody, YAP2 Antibody, Yes-associated protein 1 Antibody, Yes-associated protein 2 Antibody, YKI Antibody, YAP65 Antibody
Host	YAP1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	YAP / YAP1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	YAP / YAP1 antibody was raised against kLH-conjugated synthetic phosphopeptide (aa119-133) around Ser127 of human YAP.
Specificity:	Human YAP / YAP1
Epitope:	Ser127
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences