

ACOX1 / ACOX Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C274841 - LSBio	
CatalogID:	LS-C274841
Target:	acyl-CoA oxidase 1, palmitoyl (ACOX1)
Synonyms:	ACOX1 Antibody, ACOX Antibody, Acyl-CoA oxidase 1, palmitoyl Antibody, Acyl-coa oxidase, peroxisomal Antibody, AOX Antibody, SCOX Antibody, PALMCOX Antibody, Palmitoyl-CoA oxidase Antibody
Host	ACOX1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ACOX1 / ACOX antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ACOX1 / ACOX antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human ACOX1.
Specificity:	Human ACOX1 / ACOX.
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	