

ADK / Adenosine Kinase Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C274900 - LSBio

CatalogID:	LS-C274900
Target:	adenosine kinase (ADK)
Synonyms:	ADK Antibody, Adenosine 5-phosphotransferase Antibody, Adenosine kinase Antibody, AK Antibody
Host	ADK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ADK / Adenosine Kinase antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ADK / Adenosine Kinase antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human ADK.
Specificity:	Human ADK / Adenosine Kinase
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences