

ADRB2 Rabbit anti-Human Polyclonal (Ser364) (FITC) Antibody - LS-C274970 - LSBio

CatalogID:	LS-C274970
Target:	adrenoceptor beta 2, surface (ADRB2)
Synonyms:	ADRB2 Antibody, ADRB2R Antibody, Adrenoceptor beta 2, surface Antibody, ADRBR Antibody, B2AR Antibody, Beta-2 adrenoceptor Antibody, Beta-2 adrenergic receptor Antibody, Beta-2 adrenoreceptor Antibody, BETA2AR Antibody, BAR Antibody, Adrenergic beta-2 receptor Antibody, Catecholamine receptor Antibody
Family / Subfamily:	GPCR / Adrenergic
Host	ADRB2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ADRB2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ADRB2 antibody was raised against kLH-conjugated synthetic peptide around Ser364 of human BAR2.
Specificity:	Human ADRB2.
Epitope:	Ser364
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences