

AKT1 / AKT Rabbit anti-Human Polyclonal (Ser473) (PE) Antibody - LS-C275020 - LSBio	
<b>CatalogID:</b>	LS-C275020
<b>Target:</b>	v-akt murine thymoma viral oncogene homolog 1 (AKT1)
<b>Synonyms:</b>	AKT1 Antibody, AKT Antibody, AKT-1 Antibody, C-AKT Antibody, Pan-AKT Antibody, PKB alpha Antibody, PKBalpha Antibody, Proto-oncogene c-Akt Antibody, Protein kinase B Antibody, Rac protein kinase alpha Antibody, PRKBA Antibody, Protein kinase B alpha Antibody, RAC-PK-alpha Antibody, PKB Antibody, PKB-ALPHA Antibody, RAC Antibody, RAC-ALPHA Antibody
<b>Family / Subfamily:</b>	Protein Kinase / RAC/AKT
<b>Host</b>	AKT1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	AKT1 / AKT antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	AKT1 / AKT antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser473 of human AKT1.
<b>Specificity:</b>	Human AKT1 / AKT
<b>Epitope:</b>	Ser473
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	