

| AKT1 / AKT Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C275026 - LSBio | |
|---|---|
| CatalogID: | LS-C275026 |
| Target: | v-akt murine thymoma viral oncogene homolog 1 (AKT1) |
| Synonyms: | AKT1 Antibody, AKT Antibody, AKT-1 Antibody, C-AKT Antibody, Pan-AKT Antibody, PKB alpha Antibody, PKBalpha Antibody, Proto-oncogene c-Akt Antibody, Protein kinase B Antibody, Rac protein kinase alpha Antibody, PRKBA Antibody, Protein kinase B alpha Antibody, RAC-PK-alpha Antibody, PKB Antibody, PKB-ALPHA Antibody, RAC Antibody, RAC-ALPHA Antibody |
| Family / Subfamily: | Protein Kinase / RAC/AKT |
| Host | AKT1 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | AKT1 / AKT antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | AKT1 / AKT antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human AKT1. |
| Specificity: | Human AKT1 / AKT |
| Epitope: | C-Terminus |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/30/2014 | |
| © 2014 LifeSpan BioSciences | |