

AP2A2 Rabbit anti-Human Polyclonal (Internal) (APC) Antibody - LS-C275159 - LSBio

CatalogID:	LS-C275159
Target:	adaptor-related protein complex 2, alpha 2 subunit (AP2A2)
Synonyms:	AP2A2 Antibody, Adaptin, alpha B Antibody, ADTAB Antibody, Alpha-adaptin C Antibody, AP-2 complex subunit alpha-2 Antibody, CLAPA2 Antibody, DKFZP564D1864 Antibody, HIP-9 Antibody, HIP9 Antibody, Huntingtin yeast partner J Antibody, HYPJ Antibody, KIAA0899 Antibody, Alpha2-adaptin Antibody
Host	AP2A2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	AP2A2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP2A2 antibody was raised against kLH-conjugated synthetic peptide from internal region of human AP2A2.
Specificity:	Human AP2A2
Epitope:	Internal
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences