

ARG1 / Arginase 1 Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C275212 - LSBio

CatalogID:	LS-C275212
Target:	arginase 1 (ARG1)
Synonyms:	ARG1 Antibody, Arginase, liver Antibody, Arginase 1 Antibody, Liver-type arginase Antibody, Liver arginase Antibody, Arginase-1 Antibody, Type I arginase Antibody
Host	ARG1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ARG1 / Arginase 1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ARG1 / Arginase 1 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human ARG1.
Specificity:	Human ARG1 / Arginase 1.
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences