

BAD Rabbit anti-Human Polyclonal (Ser99) (FITC) Antibody - LS-C275405 - LSBio

CatalogID:	LS-C275405
Target:	BCL2-associated agonist of cell death (BAD)
Synonyms:	BAD Antibody, BBC2 Antibody, Bcl-2-like protein 8 Antibody, BCL2-binding component 6 Antibody, BCL2L8 Antibody, BBC6 Antibody, BCL-X/BCL-2 binding protein Antibody, Bcl-2-binding component 6 Antibody, Bcl2 antagonist of cell death Antibody, BCL2-binding protein Antibody, Bcl2-L-8 Antibody
Family / Subfamily:	Apoptosis / Not assigned-Other
Host	BAD antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	BAD antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	BAD antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser99 of human Bad.
Specificity:	Human BAD
Epitope:	Ser99
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences