

BAX Rabbit anti-Human Polyclonal (N-Terminus) (PE) Antibody - LS-C275445 - LSBio

CatalogID:	LS-C275445
Target:	BCL2-associated X protein (BAX)
Synonyms:	BAX Antibody, BCL2-associated X protein Antibody, Bcl2-L-4 Antibody, Apoptosis regulator BAX Antibody, Bcl-2-like protein 4 Antibody, BCL2L4 Antibody
Family / Subfamily:	Apoptosis / Not assigned-Other
Host	BAX antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BAX antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BAX antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region (BH3 domain) of human Bax.
Specificity:	Human Bax.
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences