

BIK Rabbit anti-Human Polyclonal (Thr33) (PE) Antibody - LS-C275477 - LSBio

CatalogID:	LS-C275477
Target:	BCL2-interacting killer (apoptosis-inducing) (BIK)
Synonyms:	BIK Antibody, Apoptosis inducer NBK Antibody, Bcl-2-interacting killer Antibody, BP4 Antibody, Apoptosis-inducing NBK Antibody, BIP1 Antibody, NBK Antibody
Host	BIK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BIK antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	BIK antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Thr33 of human Bik.
Specificity:	Human Bik
Epitope:	Thr33
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences