

ALK3 / BMPR1A Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C275544 - LSBio

<b>CatalogID:</b>	LS-C275544
<b>Target:</b>	bone morphogenetic protein receptor, type IA (BMPR1A)
<b>Synonyms:</b>	BMPR1A Antibody, 10q23del Antibody, ALK3 Antibody, ALK-3 Antibody, BMP type -1A receptor Antibody, BMPR-1A Antibody, CD292 Antibody, CD292 antigen Antibody, Activin receptor-like kinase 3 Antibody, ACVRLK3 Antibody, SKR5 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / Activin/TGFb
<b>Host</b>	BMPR1A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ALK3 / BMPR1A antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	ALK3 / BMPR1A antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human BMPR1A.
<b>Specificity:</b>	Human ALK3 / BMPR1A
<b>Epitope:</b>	N-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences