

ALK-6 / BMPR1B Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C275580 - LSBio

CatalogID:	LS-C275580
Target:	bone morphogenetic protein receptor, type IB (BMPR1B)
Synonyms:	BMPR1B Antibody, Activin receptor-like kinase 6 Antibody, ALK-6 Antibody, Bmp receptor 1b Antibody, Bmp receptor ib Antibody, BMP receptor type IB Antibody, CDw293 Antibody, CFK-43a Antibody, BMPR-1B Antibody, ALK6 Antibody, Type IB BMPR Antibody, BMP type-1B receptor Antibody, BMPR-IB Antibody, CDw293 antigen Antibody, Type IB BMP receptor Antibody
Family / Subfamily:	Protein Kinase / Activin/TGFb
Host	BMPR1B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ALK-6 / BMPR1B antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ALK-6 / BMPR1B antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human BMPR1B.
Specificity:	Human ALK-6 / BMPR1B
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences