

BRAF / B-Raf Rabbit anti-Human Polyclonal (Thr439) (FITC) Antibody - LS-C275592 - LSBio

CatalogID:	LS-C275592
Target:	v-raf murine sarcoma viral oncogene homolog B (BRAF)
Synonyms:	BRAF Antibody, 94 kDa B-raf protein Antibody, B-raf Antibody, BRAF1 Antibody, B Raf Antibody, Proto-oncogene B-Raf Antibody, p94 Antibody, B-RAF1 Antibody, NS7 Antibody, ONCOGENE BRAF Antibody, RAFB1 Antibody
Family / Subfamily:	Protein Kinase / Raf
Host	BRAF antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	BRAF / B-Raf antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	BRAF / B-Raf antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Thr439 of human BRAF.
Specificity:	Human BRAF / B-Raf
Epitope:	Thr439
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences