

BRAF / B-Raf Rabbit anti-Human Polyclonal (Internal) (FITC) Antibody - LS-C275622 - LSBio	
<b>CatalogID:</b>	LS-C275622
<b>Target:</b>	v-raf murine sarcoma viral oncogene homolog B (BRAF)
<b>Synonyms:</b>	BRAF Antibody, 94 kDa B-raf protein Antibody, B-raf Antibody, BRAF1 Antibody, B Raf Antibody, Proto-oncogene B-Raf Antibody, p94 Antibody, B-RAF1 Antibody, NS7 Antibody, ONCOGENE BRAF Antibody, RAFB1 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / Raf
<b>Host</b>	BRAF antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	BRAF / B-Raf antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	BRAF / B-Raf antibody was raised against kLH-conjugated synthetic peptide from internal region of human BRAF.
<b>Specificity:</b>	Human BRAF / B-Raf
<b>Epitope:</b>	Internal
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 9/30/2014 © 2014 LifeSpan BioSciences	