

ABL1 / c-ABL Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C275727 - LSBio	
CatalogID:	LS-C275727
Target:	c-abl oncogene 1, non-receptor tyrosine kinase (ABL1)
Synonyms:	ABL1 Antibody, Bcr/c-abl oncogene protein Antibody, Bcr-Abl Antibody, Bcr/abl Antibody, C-ABL Antibody, ABL Antibody, JTK7 Antibody, p150 Antibody, Tyrosine -protein kinase ABL1 Antibody, Proto-oncogene c-Abl Antibody, V-abl Antibody
Family / Subfamily:	Protein Kinase / Abl
Host	ABL1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ABL1 / c-ABL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ABL1 / c-ABL antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human ABL1.
Specificity:	Human ABL1 / c-ABL
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	