

**CBLC Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C275734 - LSBio**

<b>CatalogID:</b>	LS-C275734
<b>Target:</b>	Cbl proto-oncogene C, E3 ubiquitin protein ligase (CBLC)
<b>Synonyms:</b>	CBLC Antibody, CBL-3 Antibody, CBL3 Antibody, RNF57 Antibody, SH3-binding protein CBL-C Antibody, SH3-binding protein CBL-3 Antibody, CBL-c Antibody, CBL-SL Antibody, RING finger protein 57 Antibody
<b>Host</b>	CBLC antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CBLC antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CBLC antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human Cbl-c.
<b>Specificity:</b>	Human CBLC
<b>Epitope:</b>	C-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences