

CBLC Rabbit anti-Human Polyclonal (N-Terminus) (FITC) Antibody - LS-C275740 - LSBio	
CatalogID:	LS-C275740
Target:	Cbl proto-oncogene C, E3 ubiquitin protein ligase (CBLC)
Synonyms:	CBLC Antibody, CBL-3 Antibody, CBL3 Antibody, RNF57 Antibody, SH3-binding protein CBL-C Antibody, SH3-binding protein CBL-3 Antibody, CBL-c Antibody, CBL-SL Antibody, RING finger protein 57 Antibody
Host	CBLC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CBLC antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CBLC antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human Cbl-c.
Specificity:	Human CBLC
Epitope:	N-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	