

CDC25A Rabbit anti-Human Polyclonal (Ser292) (FITC) Antibody - LS-C276078 - LSBio

CatalogID:	LS-C276078
Target:	cell division cycle 25A (CDC25A)
Synonyms:	CDC25A Antibody, CDC25A2-CAG isoform Antibody, Cell division cycle 25A Antibody, CDC25A phosphatase Antibody, CDC25Hu1 Antibody, M-phase inducer phosphatase 1 Antibody, CDC25A2 Antibody
Family / Subfamily:	Protein Phosphatase / CDC25
Host	CDC25A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CDC25A antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CDC25A antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser292 of human Cdc25A.
Specificity:	Human CDC25A
Epitope:	Ser292
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences