

CHEK1 / CHK1 Rabbit anti-Human Polyclonal (Ser317) (PE) Antibody - LS-C276227 - LSBio

CatalogID:	LS-C276227
Target:	checkpoint kinase 1 (CHEK1)
Synonyms:	CHEK1 Antibody, CHK1 checkpoint homolog Antibody, Cell cycle checkpoint kinase Antibody, Checkpoint kinase 1 Antibody, CHK1 Antibody, Chk1-S Antibody, Protein kinase chk1 Antibody, Checkpoint kinase-1 Antibody
Family / Subfamily:	Protein Kinase / NIM1
Host	CHEK1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CHEK1 / CHK1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CHEK1 / CHK1 antibody was raised against kLH-conjugated synthetic phosphopeptide from amino acids around Ser317 of human CHK1.
Specificity:	Human CHEK1 / CHK1
Epitope:	Ser317
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences