

CLTA / LCA Rabbit anti-Human Polyclonal (Internal) (Biotin) Antibody - LS-C276289 - LSBio	
CatalogID:	LS-C276289
Target:	clathrin, light chain A (CLTA)
Synonyms:	CLTA Antibody, Clathrin light chain A Antibody, Clathrin, light chain A Antibody, LCA Antibody
Host	CLTA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Biotin
Immunogen Species:	CLTA / LCA antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CLTA / LCA antibody was raised against kLH-conjugated synthetic peptide from internal region of human CLTA.
Specificity:	Human CLTA / LCA
Epitope:	Internal
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	